



STORAGE INSTRUCTIONS: TARGA® II HF CARTRIDGES

A. SHELF LIFE

The cartridges leave the factory stored in a glycerin and water solution. From the date of shipment, the cartridge shelf life in the original packaging is one year. This assumes that the sealed bags have not been opened or damaged, and the conditions shown in Sections B&C are met. Consult KMS if cartridges have been stored for a longer period of time.

B. HANDLING

1. The boxes should not be dropped or crushed.
2. The cartridges should be stored horizontally.

C. CARTRIDGE STORAGE CONDITIONS

1. Keep the cartridges in their original packaging and store in a horizontal position.
2. Store indoors out of direct sunlight.
3. Do not allow the cartridges to freeze, and they should be kept at less than 40°C (104°F).
4. The preferred long-term storage temperature range is 10°C to 30°C (50 to 85°F), at relative humidity below 70%.

D. USED CARTRIDGE STORAGE SOLUTION

1. **Short-term shut-down (up to 3 days):** Cartridges should be thoroughly cleaned and left in water at 5-25°C (41 to 77°F).
2. **Long-term shut-down (3 days or longer):** cartridges should be cleaned, rinsed and impregnated before storage with one of the following solutions:
 - a. Glycerin @ 80 – 100% (best)
 - b. Phosphoric acid @ pH 2 – 3
 - c. Benzoic acid or sodium benzoate @ 1,000 ppm
 - d. Sodium bisulfite or sodium metabisulfite @ 1,000 – 5,000 ppm

Cartridges stored in sodium bisulfite or sodium metabisulfite should be flushed with clean water and impregnated with fresh solution at six-month intervals. Once impregnated, the cartridges should be left on the system with all valves closed or removed from the system and sealed in plastic bags. Storage conditions described above for new cartridges should also be used for storage of used cartridges.

E. POST-STORAGE CLEANING PROCEDURE

All cartridges should be drained, rinsed and cleaned according to KMS' recommended pre-startup procedures prior to restarting the system on production.

Koch Membrane Systems, Inc., www.kochmembrane.com

Corporate Headquarters: 850 Main Street, Wilmington, Massachusetts 01887-3388, US, Tel. Toll Free: 1-888-677-5624, Telephone: 1-978-694-7000, Fax: 1-978-657-5208

European Headquarters: Koch Chemical Technology Group Ltd., Units 3-6, Frank Foley Way, Stafford ST16 2ST, GB, Telephone: +44-178-527-2500, Fax: +44-178-522-3149

• Aachen DE • Lyon FR • Madrid ES • Milan IT • Wijnegem BE • Beijing & Shanghai CN • Mumbai & New Delhi IN • Melbourne & Sydney AU • Singapore • Sao Paulo BR • Dubai UAE •

For related trademark and patent information, visit www.kochmembrane.com/legal

Koch Membrane Systems, Inc. is a Koch Chemical Technology Group, LLC company.

© 2015 Koch Membrane Systems, Inc. All rights reserved worldwide.

01/15 Rev 15-1